



UNITE DE SOINS
ET DE **CARDIOLOGIE**
INTERVENTIONNELLE



NOUVELLES CLINIQUES NANTAISES



Recent Advances in Electrical Therapies



Nice, Oct 2013

Gras, MD, Nantes, France

Recent Advances in Electrical Therapies

- **Leadless PM: LLP**
- Subcutaneous Lead ICD: S-ICD
- Recent CRT Recommendations
- Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP
- LV Endocardial Pacing: LVEP
- Vagal Nerve stimulation in HF: VNS
- Left Atrial Appendage Occlusion: LAAO

Unmet Needs in Cardiac Pacing & ICD


- **Pacing Lead failure**

- *Up to 21% within 10 years after PM implantation*




- **ICD Lead failure**

- *Lead replacement is mandatory in 38% within 8 years after ICD implantation*



- Lead Fracture & Dislodgement
- Vasculatures and Cardiac Infections
- PNO Hemothorax, Cardiac Perforation
- Tricuspid Valve Injury
- Vascular Thrombosis



Subclavian Vein/IT Occlusion



PM Lead Issues in Clinical Practice



Special Article

Totally Self-Contained Intracardiac Pacemaker*

J. WILLIAM SPICKLER, PH.D., NED S. RASOR, PH.D., PAUL KEZDI, M.D.,
S. N. MISRA, M.D., K. E. ROBINS, P.E., AND CHARLES LeBOEUF, P.E.

SUMMARY

Recent developments in miniature long-life power sources and electronics, such as nuclear batteries and integrated circuits make feasible a new generation of pacemakers, the intracardiac pacemaker (IC), i.e., a completely self-contained pacemaker implanted inside the right ventricle by transvenous insertion. Since the IC pacemaker eliminates all leads, problems associated with the leads such

circuits have been improved substantially. In addition, the development of the endocatheter electrode has broadened the range of operative procedures to include a portion of the patient population. Two major problems that still exist with conventional pacemakers are perforation or dislocation of the transvenous electrode and the short life of the batteries that are presently used. In addition, there is a certain physical and psychological discomfort involved with the




Fig. 2. Some early unsatisfactory dummy capsules




Fig. 4. Intracardiac pacemaker with catheter for transvenous insertion.





Fig. 8. Nuclear-powered intracardiac pacemaker.



Leadless Pacemaker Potential Benefits and Risks



Miniaturized, Leadless
VVIR Pacer



Steerable Sheath/Catheter



Reduced Invasiveness

- *No surgery, Femoral venous access*
- *Less radiation exposure (femoral)*
- *More cosmetic for patient (“invisible”)*

Improved Efficiency

- *Fewer complications (no lead or subQ device)*
- *MRI conditional*

More Cost-Effective Therapy

- *Reduced length of hospital stay*
- *Fewer acute and chronic complications*

Risks (to be evaluated)

- *Device dislodgement, Migration*
- *AV Dyssynchrony*
- *Explantation ??*

Leadless LV Pacing

- Upgrade existing implanted systems
 - Works with any PM or ICD
- Simple co-implant
 - Transvenous right side system
 - Wireless left side system








Figure 2 An example of clinically determined acoustic windows in 4 body positions (in red with the patient lying supine; in green with 30° right tilt; in yellow with 30° left tilt; in purple with 30° upright tilt) superimposed on the CT-determined acoustic window (in light blue with the patient lying supine and during end inspiration) on 3D reconstruction CT of the thorax. 3D = three-dimensional; CT = computed tomography.

Recent Advances in Electrical Therapies

- Leadless PM: LLP
- **Subcutaneous Lead ICD: S-ICD**
- Recent CRT Recommendations
- Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP
- LV Endocardial Pacing: LVEP
- Vagal Nerve stimulation in HF: VNS
- Left Atrial Appendage Occlusion: LAAO

“Minimally Invasive ICD” – Pediatric Case Reports: Children with congenital heart disease





- SQ array around left thorax
- Abdominal active can generator
- Epicardial pace-sense lead



DFT \leq 14 - 20 joules

Unique challenges of ICD implant
in children due to size, growth,
& frequent congenital heart disease
with increased risk paradoxical embolism

Stephenson J CV EP 2006;17:41-6
Berul PACE 2001;24:1789-4;
Gradaus JCE 2001;12:356-60




N Engl J Med 2010;363:36-44.






S-ICD Implantation Procedure

- . Tripolar parasternal electrode is positioned parallel to and 1 to 2 cm to the left of the sternal midline
- . Distal sensing electrode placed close to the manubriosternal junction
- . Proximal sensing electrode placed adjacent to the xiphoid process
- . PG is positioned over the sixth rib between the midaxillary line and the anterior axillary line.
- . No fluoroscopy is required



The electrode has an 8-cm shocking coil, flanked by 2 sensing electrodes.



S-ICD Implantation Procedure







Sensing the Subcutaneous Signal

Automatic Setup Analyzes:

- QRS and T-wave amplitudes
- QRS:T-wave ratios (amplitude and timing)
- Patient posture
- QRS morphology and width



Screening prior to and post S-ICD implantation



Variant Patient Posture


What About Acceptance of S-ICD?

- Inappropriate shocks due to oversensing (80 J output)
- Potential erosion/migration issues (large device 145 g; 58 cc)
- Tunnelling along rib margin and parasternally
- No painless therapies (ATP)
- No bradycardia support
 - ~5-10% of prophylactic ICD pts will develop concomitant bradycardia
 - Aggressive use of BB could exacerbate this condition (HF population)
- No heart failure monitors
- No remote patient follow up


INSIGHT™ Rhythm Discrimination

3 methods to correctly identify & classify the S-ECG rhythm:

1. Static Morphology Analysis
2. Dynamic Morphology Analysis
3. QRS Width Analysis




INSIGHT™ Rhythm Discrimination: PVT/VF



1. Static morphology (red line) and
2. Dynamic morphology (green line) fall during PVT/VF.
3. QRS Width (pink) increases during Vent tachyarrhythmias

INSIGHT™ Rhythm Discrimination: VT



1. Static morphology (red line) drops during VT while
2. Dynamic morphology (green line) remains high due to beat-to-beat similarity.

S-ICD Clinical Cases

Device Settings

Therapy: ON

Shock Zone: 240 bpm

Conditional Shock Zone: 220 bpm

Post Shock Pacing: ON

Gain Setting: 1X

Sensing Configuration: Alternate

S = Sense

P = Pace

N = Noise

T = Tachy Detection



C = Charge Start

. = Discard

⚡ = Shock

▬ = Episode End

TREATED EPISODE 003: 05/26/2011
SHOCK IMPEDANCE = 48 Ohms



Inappropriate shocks due to T wave oversensing

The Entirely Subcutaneous Implantable Cardioverter-Defibrillator


118 patients, FU 18 Months

8 pts: 45 successful appropriate shocks

15 pts (13%) inappropriate shocks:

T Wave oversensing, myopotentials, double counting, AF


16 (14%) Complications:



Recent Advances in Electrical Therapies

- Leadless PM: LLP
- Subcutaneous Lead ICD: S-ICD
- **Recent CRT Recommendations**
- Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP
- LV Endocardial Pacing: LVEP
- Vagal Nerve stimulation in HF: VNS
- Left Atrial Appendage Occlusion: LAAO

CRT-D - Units per million inhabitants



Source population data: OECD

Units - Eucomed based on reports from major manufacturers

* Europe represents total of listed countries (N/A countries excluded)



2013 ESC Guidelines on cardiac pacing and cardiac resynchronization therapy

The Task Force on cardiac pacing and resynchronization therapy of the European Society of Cardiology (ESC). Developed in collaboration with the European Heart Rhythm Association (EHRA).

Authors/Task Force Members: Michele Brignole (Chairperson) (Italy)*, Angelo Auricchio (Switzerland), Gonzalo Baron-Esquivias (Spain), Pierre Bordachar (France), Giuseppe Boriani (Italy), Ole-A Breithardt (Germany), John Cleland (UK), Jean-Claude Deharo (France), Victoria Delgado (Netherlands), Perry M. Elliott (UK), Bulent Gorenek (Turkey), Carsten W. Israel (Germany), Christophe Leclercq (France), Cecilia Linde (Sweden), Lluís Mont (Spain), Luigi Padeletti (Italy), Richard Sutton (UK), Panos E. Vardas (Greece)

ESC Committee for Practice Guidelines (CPG): Jose Luis Zamorano (Chairperson) (Spain), Stephan Achenbach (Germany), Helmut Baumgartner (Germany), Jeroen J. Bax (Netherlands), Héctor Bueno (Spain), Veronica Dean (France), Christi Deaton (UK), Cetin Erol (Turkey), Robert Fagard (Belgium), Roberto Ferrari (Italy), David Hasdai (Israel), Arno W. Hoes (Netherlands), Paulus Kirchhof (Germany/UK), Juhani Knuuti (Finland), Philippe Kolh (Belgium), Patrizio Lancellotti (Belgium), Ales Linhart (Czech Republic), Petros Nihoyannopoulos (UK), Massimo F. Piepoli (Italy), Piotr Ponikowski (Poland), Per Anton Sirnes (Norway), Juan Luis Tamargo (Spain), Michal Tendera (Poland), Adam Torbicki (Poland), William Wijns (Belgium), Stephan Windecker (Switzerland).

Predictive value of QRS duration in NYHA class I-II patients

MADIT CRT

QRS duration


<150 msec

147/645

P for interaction= 0.001

≥150 msec

225/1175



REVERSE 24-months

QRS interval

<152 ms


120

0.80 (0.35-1.85)

≥152 ms

142

0.23 (0.10-0.56)



RAFT

QRS duration


Intrinsic QRS <150 msec

248/627

Intrinsic QRS ≥150 msec

359/1036

P for interaction= 0.003



Effectiveness of Cardiac Resynchronization Therapy by QRS Morphology in the Multicenter Automatic Defibrillator Implantation Trial Cardiac Resynchronization Therapy (MADIT-CRT)


Wojciech Zareba, Helmut Klein, Iwona Cygankiewicz, W. Jackson Hall, Scott McNitt, Mary Brown, David Cannom, James P. Daubert, Michael Eldar, Michael R. Gold, Jeffrey J. Goldberger, Ilan Goldenberg, Edgar Lichstein, Heinz Pitschner, Mayer Rashtian, Scott Solomon, Sami Viskin, Paul Wang, Arthur J. Moss and on behalf of the MADIT-CRT Investigators

Circulation 2011;123:1061-1072; originally published online Feb 28, 2011;


MADIT CRT *Primary endpoint: All-cause death or Heart Failure event*

Predictive value of QRS morphology in NYHA class I-II pts


LBBB



Non-LBBB





REVERSE 5-years: Mortality by QRS Morphology




Patients in SR + LBBB

NYHA-II, III et IV (ambulatory),
OMT, EF ≤ 35%



Patients in SR + no LBBB

NYHA-II, III et IV (ambulatory),
OMT, EF \leq 35%




Cardiac-Resynchronization Therapy in Heart Failure with a Narrow QRS Complex

Frank Ruschitzka, M.D., William T. Abraham, M.D., Jagmeet P. Singh, M.D., Ph.D., Jeroen J. Bax, M.D., Ph.D., Jeffrey S. Borer, M.D., Josep Brugada, M.D., Ph.D., Kenneth Dickstein, M.D., Ph.D., Ian Ford, M.D., Ph.D., John Gorcsan III, M.D., Daniel Gras, M.D., Henry Krum, M.B., B.S., Ph.D., Peter Sogaard, M.D., D.M.Sc., and Johannes Holzmeister, M.D., for the EchoCRT Study Group*

The NEW ENGLAND JOURNAL of MEDICINE


A Primary Composite Outcome



No. at Risk

CRT	404	297	223	155	103	65	42	19
Control	405	302	236	166	119	71	44	15

B Death from Any Cause



No. at Risk

CRT	404	334	267	199	132	84	56	25
Control	405	335	269	195	141	87	62	27

2013 ESC Guidelines on cardiac pacing and cardiac resynchronization therapy




Figure 10 Indication for atrioventricular junction (AVJ) ablation in patients with symptomatic permanent atrial fibrillation (AF) and optimal pharmacological therapy. BiV = biventricular; CRT = cardiac resynchronization therapy; EF = ejection fraction; HR = heart rate; ICD = implantable cardioverter defibrillator; NYHA = New York Heart Association.

HF Patients with conventional indication for Cardiac Pacing

Upgrade from conventional PM or ICD

De Novo CRT

Indicated in case of

- EF < 35%, OMT
- classe NYHA-III- IV (ambulatory)
- High % of VP

Classe I, Level B

Indicated in case of:

- Reduced EF (<50%)
- Expected high % of VP

Classe IIa, Level B

Need for Pacing

NYHA III/IV

NYHA II

EF ≤ 35%

EF ≤ 35%

any QRS

any QRS



Class IIa C

CRT-P
CRT-D

Class IIb C

CRT-P
CRT-D

Choice of pacing mode (and cardiac resynchronization therapy optimization)

Recommendations	Class ^a	Level ^b	Ref. ^c
1) The goal of CRT should be to achieve BiV pacing as close to 100% as possible since the survival benefit and reduction in hospitalization are strongly associated with an increasing percentage of BiV pacing.	IIa	B	67–69
2) Apical position of the LV lead should be avoided when possible.	IIa	B	70–72
3) LV lead placement may be targeted at the latest activated LV segment.	IIb	B	73

Left Ventricular Lead Position and Clinical Outcome in the Multicenter Automatic Defibrillator Implantation Trial–Cardiac Resynchronization Therapy (MADIT-CRT) Trial

Jagmeet P. Singh, MD, DPhil*; Helmut U. Klein, MD*; David T. Huang, MD; Sven Reek, MD; Malte Kuniss, MD; Aurelio Quesada, MD; Alon Barsheshet, MD; David Cannom, MD; Ilan Goldenberg, MD; Scott McNitt, MS; James P. Daubert, MD; Wojciech Zareba, MD; Arthur J. Moss, MD






Figure 3: Percent change in reverse LV remodeling parameters at 6-month follow-up: LVEF (A), LV EDD (B), and LV ESD (C). Data are presented as mean and standard error of the mean.

The Target Study



All Cause Mortality following CRT in the TARGET and Control Groups



Combined Endpoint of Death and Heart Failure Related Hospitalisation between the TARGET and Control Groups

Echocardiography-Guided Left Ventricular Lead Placement for Cardiac Resynchronization Therapy: Results of the Speckle Tracking Assisted Resynchronization Therapy for Electrode Region (STARTER) Trial

Samir Saba MD, Josef Marek MD, David Schwartzman MD, Sandeep Jain MD, Evan Adelstein MD, Pamela White RN, Olusegun A. Oyenuga MD, Tetsuari Onishi MD, Prem Soman MD, John Gorcsan III MD University of Pittsburgh Medical Center, Pittsburgh, PA

Address correspondence to: Samir Saba, MD University of Pittsburgh Medical Center 200 Lothrop Street, Suite B-535 Pittsburgh, PA 15213-2582 sabas@upmc.edu

In Press: Circulation Heart Failure 2013

Recent Advances in Electrical Therapies



- Leadless PM: LLP
- Subcutaneous Lead ICD: S-ICD
- Recent CRT Recommendations
- **Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP**
- LV Endocardial Pacing: LVEP
- Vagal Nerve stimulation in HF: VNS
- Left Atrial Appendage Occlusion: LAAO

Atrial based, Dual Site LV, RV Pacing


Clinical Trials

Addition of a Second LV Pacing Site in CRT Nonresponders Rationale and Design of the Multicenter Randomized V³ Trial



PIERRE BORDACHAR, MD,¹ CHRISTINE ALONSO, MD,² FREDERIC ANSELME, MD,³ SERGE BOVEDA, MD,⁴
PASCAL DEFAYE, MD,⁵ STEPHANE GARRIGUE, MD,⁶ DANIEL GRAS, MD,⁷ DIDIER KLUG, MD,⁸ OLIVIER PIOT, MD,⁹
NICOLAS SADOUL, MD,¹⁰ AND CHRISTOPHE LECLERCQ, MD¹¹




Atrial based, Dual Site LV, RV Pacing (ongoing V3 Trial)






Interest of Additional LV Lead during CRT





Right Subclavian vein Occlusion DDD-PM Upgrade







Additional Hemodynamic Impact of MPP


Subselection of Lateral Cardiac Vein during Quadripolar LV Lead Implant



Quadripolar LV Lead in case of LSVC



Goals of Quadripolar vs Bipolar LV Pacing



- Offering Multiple Pacing Vectors to individually adjust LV Pacing in Long-Term
 - Less Phrenic Nerve Stimulation
 - Optimal Pacing Thresholds (multiple vectors)
 - Pacing Vectors to Optimize Hemodynamics
 - Better Lead Stability, Less Lead revision
 - Similar LV Lead Implantation





Recent Advances in Electrical Therapies

- Leadless PM: LLP
- Subcutaneous Lead ICD: S-ICD
- Recent CRT Recommendations
- Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP
- **LV Endocardial Pacing: LVEP**
- Vagal Nerve stimulation in HF: VNS
- Left Atrial Appendage Occlusion: LAAO

Endocardial vs Epicardial CRT provides:

Better LV Filling and Systolic Performance

More Homogenous Resynchronization



Prinzen et al

Van Deursen, Circ Arrhythmia
Electrophysiol. 2009;2:580-587

LV Endocardial or Triventricular Pacing to Optimize CRT in a Chronic Canine Model of Ischemic HF

Bordachar Am J Physiol Heart Circ Physiol 2012.

Table 1. Comparison of Hemodynamics during LV and BiV pacing (n=8)


Table 3. Distribution of optimal sites during epicardial and endocardial pacing

Site \ Pacing	LV EPI	LV ENDO	BiV EPI	BiV ENDO
LVA	0%	38%	12%	12%
LVL	50%	12%	13%	25%
LVB	25%	38%	50%	63%
LVPOST	25%	12%	25%	0%

LVA: LV apex; LVL: LV lateral; LVB: LV basal; LVPOST: LV posterior

Optimal Left Ventricular Endocardial Pacing Sites for Cardiac Resynchronization Therapy in Patients With Ischemic Cardiomyopathy


David D. Spragg, MD,* Jun Dong, MD, PHD,*† Barry J. Fetters, MS,* Robert Helm, MD,* Joseph E. Marine, MD,* Alan Cheng, MD,* Charles A. Henrikson, MD,* David A. Kass, MD,* Ronald D. Berger, MD, PHD*




CRT delivered at best LV endocardial sites is more effective than via pre-implanted coronary sinus lead pacing. The location of optimal LV endocardial pacing varies among patients with ICM, and individual tailoring may improve CRT efficacy in such patients. (J Am Coll Cardiol 2010;56:774-81) © 2010 by the American College of

TEE Evaluation before Transseptal Puncture

Location of Fossa Ovalis



VE IM 0,5 26/04/2011
04-26-110414 N.C.N. CARDIOLOGIE ITm 0,0 11:12:49



Ongoing Evaluation

LV Endo Pacing in Non CRT Responder

LV Lead placement
under TEE guidance


PHILIPS HERVE IM 0,5 26/04/2
11-04-26-110414 N.C.N. CARDIOLOGIE ITm 0,0 13:09:5

ETO
X7-2t
44Hz
14cm

T. Pat.: 37,0 °C
T ETO : 38,3 °C
0 84 180



2D
Gén
Gn 70
C 50
4 / 2 / 0
75 mm/s

G
P R
3,0 8,0



Ongoing Evaluation

LV Endo Pacing in Non CRT Responder




AP View


Ongoing Evaluation

LAO View

LV endocardial Pacing during CRT





AP View



LAO View

PNS still happens during LV endo Pacing !!



Echo prior to Transseptal LV Lead implant

FE VG (sim.A4C)

19 %

Vol. Eject.(simpson A4C)

25 ml

2 VGs Long(A4C)

8.5 cm

Vol.VG.tS(sim.4cav)



106 ml

1 VGd Long(4cav)


9.2 cm

Vol.VG.tD(sim.4cav)


131 ml



08/04/2011 12:01:41




LV EF Echo Evaluation (Simson)



Permanent LV Endocardial Pacing in Clinical Practice

- Expected Benefits of LVEP:
 - Local Recruitment, V propagation
 - Optimized LV Pacing Location
 - Lower risk of PNS, Better PT
- Potential Side Effects:
 - TE events
 - Impact on MR
 - Lead Extraction, arrhythmias ..
- Indications:
 - Non Responders
 - Failure to Classical Approach
 - First Line option ?



- Next Steps
 - Clinical Studies
 - Improved Implant Tools
 - New LVEP Lead Design

Recent Advances in Electrical Therapies

- Leadless PM: LLP
- Subcutaneous Lead ICD: S-ICD
- Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP
- LV Endocardial Pacing: LVEP
- **Vagal Nerve stimulation in HF: VNS**

Downloaded from heart.bmj.com on January 10, 2013 - Published by group.bmj.com

Vagal stimulation for the treatment of heart failure: a translational success story

Peter J Schwartz^{1,2,3,4}

inhibition of cardiac activity (ie, low the heart from sympathetic nervous system home message stimulates—as studies—the in also and significant synergistic effects a significant wa


Device-based Neuromodulation Therapy for HF

Current Investigational Approaches

SCS

VNS

BRS



Spinal Cord Stimulation

(SCS): SCS generator implant in abdomen or paraspinal region with a lead placed in dorsal epidural space between T1-T4


Vagal Nerve Stimulation

(VNS): VN stimulator placed in right subpectoral region with a lead placed in RV and VN stimulating lead tunneled to cervical vagus region


Baroreflex Stimulation

(BRS): Baroreflex stimulation generator placed in right subpectoral region with bilateral stimulation leads tunneled to the carotid baroreceptor region


Impact of Vagal Stimulation



Central Effects




Cardiac Effects



Stimulating the cervical vagus nerve activates regions in the brain and in the heart

Vagal stimulation is indicated for epilepsy.




Please note:

Epilepsy treatment uses mainly left side vagus stimulation

Cyberonics has implanted >60,000 devices. Safety profile established.



Heart failure as an autonomic nervous system dysfunction

Takuya Kishi (MD, PhD)*



Principal influence on sympathetic outflow


VNS in HF Canine High Rate Pacing Model



Adapted from Zhang Y, Circ Heart Fail 2009; 2:692-699

VNS in HF Canine High Rate Pacing Model



	Baseline		4-wk Pacing		8-wk Pacing	
	Control	VNS	Control	VNS	Control	VNS
RR, ms	510.7±77.0	514.5±61.5	394.8±36.7	428.8±55.7	407.1±47.2	451.0±76.1
SDNN, ms	84.2±21.7	86.6±21.8	23.2±5.9	36.6±5.1*	28.7±8.0	42.2±7.4*
RMSSD, ms	68.9±10.6	69.3±17.2	17.0±4.6	31.0±6.1*	22.1±5.3	37.2±7.1*
LF, norm	35.2±12.5	36.2±12.3	72.1±8.6	55.6±6.1*	65.3±10.3	53.2±9.6*
HF, norm	64.8±12.5	63.8±12.3	27.9±8.6	44.4±6.1*	34.7±10.3	46.8±9.6*
LF/HF	0.70±0.33	0.63±0.34	3.03±1.79	1.29±0.33*	2.23±1.46	1.22±0.75



Adapted from Zhang Y, Circ Heart Fail 2009; 2:692-699

Chronic vagus nerve stimulation: a new and promising therapeutic approach for chronic heart failure

Gaetano M. De Ferrari^{1*}, Harry J.G.M. Crijns², Martin Borggrefe³, Goran Milasinovic⁴, Jan Smid⁵, Markus Zabel⁶, Antonello Gavazzi⁷, Antonio Sanzo¹, Robert Dennert³, Juergen Kuschyk⁴, Srdjan Raspopovic⁵, Helmut Klein^{6,8}, Karl Swedberg⁹, and Peter J. Schwartz^{1,10,11,12,13}, for the CardioFit Multicenter Trial Investigators






Figure 1 New York Heart Association class. Individual New York Heart Association class behaviour of each patient throughout the study. Red broken lines correspond to dead patients. Each patient is characterized by a single line.

Data limited to the 21 patients of the second phase of the study show a behaviour identical to the total cohort, but a non-significant decrease in LV systolic volume (see Supplementary material online, Table S1). Comparing baseline with 1-year evaluation in the 23 patients who completed this follow-up revealed maintenance and even magnification of the favourable effects of vagal stimulation (especially LVEF from 21 to 34%, Table 5).

Heart rate variability tended to increase during the study and the change in pNN50 was statistically significant (Table 4).




Figure 2 Echocardiographic evaluation. Left ventricular volume indexes (A) and ejection fraction (B) throughout the study. LVEDVI: left ventricular end-diastolic volume index. LVESVI: left ventricular end-systolic volume index.

CardioFit Multicenter Trial - Safety

- Non SAEs related to the device: n = 19

Pain at stimulation site: 6

Cough: 5

Dysphonia: 4

Mandibular pain: 3

Stimulus artefact on ECG: 1

- Non SAEs not related to the device: n = 11

NECTAR-HF Study: Protocol Overview

- **Study Design**


- Single-blind, placebo controlled, randomized 2:1 (therapy/control)
- Multicentre (European sites)
- Control patients crossed over to therapy at 6M follow-up & followed for safety through 18 months

- **Sample Size**

- 250 pts screened for eligibility
- 96 pts implanted with the system

- **Patient Population:**

- NYHA class III HF pts
- Ejection fraction of $\leq 35\%$
- Not CRT candidate, QRS ≤ 130 ms



Recent Advances in Electrical Therapies

- Leadless PM: LLP
- Subcutaneous Lead ICD: S-ICD
- Recent CRT Recommendations
- Multisite, Quadripolar LV Pacing: MSP, 4P-LVP, MPP
- LV Endocardial Pacing: LVEP
- Vagal Nerve stimulation in HF: VNS
- **Left Atrial Appendage Occlusion: LAAO**

NEW ENTRIES IN THE ESC GUIDELINES FOR THE MANAGEMENT OF ATRIAL FIBRILLATION (AF)¹

Recommendation for Left Atrial Appendage Occlusion (LAAO)

The European Society of Cardiology has released the 2012 focused update on the guidelines for the management of atrial fibrillation. This update includes for the first time a recommendation for the use of left atrial appendage occlusion.

The update was developed with contribution from the European Heart Rhythm Association (EHRA).

Recommendations for LAA closure/occlusion/excision

Recommendations	Class ^a	Level ^b	Ref ^c
Interventional, percutaneous LAA closure may be considered in patients with a high stroke risk and contraindications for long-term oral anticoagulation.	IIb	B	115, 118
Surgical excision of the LAA may be considered in patients undergoing open heart surgery.	IIb	C	

LAA = left atrial appendage.




^aClass of recommendation.



^bLevel of evidence.

^cReferences.

- "Interventional percutaneous occlusion/closure of the LAA has a role in patients with thromboembolic risk who cannot be managed in the long-term using any form of OAC."¹

Source: Camm A, Lip G, De Caterina R, et al. "2012 focused update of the ESC Guidelines for the management of atrial fibrillation." *Eur Heart J*. Presented ahead of publication European Society of Cardiology Congress; 2012 August 25-29; Munich, Germany.











PHILIPS PICARD
13-02-04-115754

IM 0,5 04/02/2013 RD
ITm 0,0 12:41:47 :04-115754

IM 0,5 04/02/2013
ITm 0,9 12:48:18







Recent Advances in Electrical Therapies

1. Quadripolar LV Pacing: PNS, PT, MPP
2. Left Atrial Appendage Occlusion
3. Leadless PM, Subcutaneous Lead ICD
4. HF Refractory to Conventional XX
 - . LV Endocardial Pacing
 - . Multisite Ventricular Pacing
 - . Vagal Nerve stimulation in HF